

NEXT GENERATION BIOCERAMIC SEALER FROM THE INVENTOR OF WHITE PROROOT®* MTA1

-000mr Sther

Cortest and the second second

ONE SEALER TO STREAMLINE YOUR PRACTICE

NOT ALL CASES CAN BE TREATED THE SAME WAY.

NeoSEALER Flo is optimized for both warm vertical and single cone obturation techniques for maximum efficiency and savings.

FORMULATED WITH **MORE BIOACTIVE CEMENT THAN ORIGINAL BIOCERAMIC SEALER²**

NeoSEALER Flo is thoughtfully designed to promote setting and maximize bioactivity.

MTA-BASED BIOACTIVITY

Hydroxyapatite forms to ensure bioactive sealing.

DIMENSIONALLY STABLE

Unlike resin-based sealers, NeoSEALER Flo is dimensionally stable.

ADA 57 std. <1% shrinkage & <0.1% expansion

I love the way Avalon Biomed NeoSEALER Flo handles. It coats the canal walls consistently and clean up in the chamber after obturation is easy.

> - Ryan Walsh, DDS, MS Endodontist


RESIN-FREE

For maximum bioactivity.

NON-STAINING

- ✓ Won't discolor teeth EVER.
- ✓ NaOCI won't cause discoloration.

EXTENDED WORKING TIME



** Sets in vivo in the presence of moisture

HIGH FLOW

Fine particle size allows excellent flow to fill and seal lateral canals.

HIGH RADIOPACITY





Courtesy of Alejandro Podolsky, DDS Courtesy of Ron Hill, DDS, MSD Endodontist

Endodontist

* PROROOT®, DENTSPLY® and DENTSPLY SIRONA® are registered trademarks of DENTSPLY SIRONA Inc 1 Patented under US Patent No. 7,892,342 owned by DENTSPLY International Inc.

2 NeoSEALER Flo contains 6% more bioactive cement by weight compared to EndoSequence BC Sealer (in-house data obtained through XRD analysis) 3 | Dent. 2014 May; 42(5):517-33. A review of the bioactivity of hydraulic calcium silicate cements. Niu L.N., Jiao K, Wang TD, Zhang W, Camilleri J, Bergeron BE, Feng HL, Mao J, Chen JH, Pashley DH, Tay FR.


HIGH pH

To promote osteogenic response.

MORE VALUE BUILT-IN, AT LESS COST

Our unique Flex Flo Tip[™] design combines a flexible cannula with an integrated stopper attached to a minimal waste barrel for highly accurate placement.

Flex Flo Tip[™] More Precision. 81% Less Waste.




These Flex Flo Tips are amazing! The combination of the flexible cannula and the rubber stopper increases user control allowing accurate and gentle guided placement of the bioceramic sealer into the canal. The tips make endodontics much more predictable and accurate.


> - Brian Minavi, DDS, MS Endodontist

UNLIKE INERT RESIN AND ZOE MATERIALS, MTA PROVIDES STATE OF THE ART TECHNOLOGY TO PROMOTE BIOACTIVITY AND SUPPORT HEALING

- After placement of MTA in the tooth, body fluids infiltrate the paste to cause setting.
- During this process calcium and hydroxide ions are released from the MTA, promoting hydroxyapatite (HA) formation on the MTA surfaces.
- When HA coats the surface, it hides the underlying MTA to minimize foreign body reactions and support healing responses of the pulpal or periapical tissue.



Calcium phosphate precipitated crystals scattered on the hydrated tricalcium silicate cement surface³.



High magnification image of the calcium phosphate crvstals³.

NEXT GENERATION BIOCERAMIC SEALER

For years, clinicians have settled for using inert resin and ZOE sealers or paid too much for premixed bioceramic sealers. **Not anymore.**

At Avalon Biomed, we believe innovation should improve patients' lives and reduce health care costs. That's why we designed NeoSEALER Flo to be efficient, effective and the most economical premixed bioceramic sealer on the market today.

From the beginning, our goal has been to save more teeth for more people. That takes constant improvement in technology, without resting on past success.

> – Carolyn Primus, PhD Inventor of White ProRoot®[®] MTA', Grey MTA Plus®, NeoMTA Plus®, NeoMTA® 2, NeoSEALER™ Flo, NeoPUTTY™



ARMAMENTARIUM with NeoPUTTY, the most radiopaque bioactive

bioceramic putty in its class.

NeoPUTTY is a resin-free bioactive bioceramic cement with superior handling properties. Its firm, low-tack consistency, wash-out resistance and bioactivity make it the premier putty for 12 of your most important procedures.

AVAILABLE IN:

Item Code	Item Description
ANSSK	NeoSEALER Flo - 2.2 gm Starter Kit (with 20 Flex Flo Tips)
ANSPK	NeoSEALER Flo - 4.4 gm Professional Kit (with 40 Flex Flo Tips)
ANSPCK	Combo Kit - NeoSEALER Flo 2.2 gm Starter Kit + NeoPUTTY 1.2 gm Professional Kit (with 20 Flex Flo Tips)
FFT20	Flex Flo Tips (20/pack)
FFT40	Flex Flo Tips (40/pack)

PLACE YOUR ORDER

and see what's possible with the world's most versatile bioactive bioceramics at

AVALONBIOMED.COM

Or email us to find out more at

INFO@AVALONBIOMED.COM



MKT20016